



## Microcontact Printing of Surface Proteins

**Prepared by:** Evren U. Azeloglu

**Version:** 1.0

**Date:** 01-Dec-2014

### Materials

- PDMS, Polydimethylsiloxane (Slygard 184, Fisher, Cat: NC9644388)
- Glass coverslips, such as 22x22-1 (Fisher, Cat: 12-548-B)
- Style 2-A Forceps (EMS, Cat: 72919-2A)
- Tween-20 (Sigma, Cat: P1379)
- Pluronic Solution (0.2% in diH<sub>2</sub>O diluted from 1% stock; prepared ahead of time)
  - △ Pluronic F-127 (Sigma, Cat: P2443)
- Liquid Silanization Solution (5% silane in benzene)
  - △ 1,2 dichlorobenzene, anhydrous (Sigma, Cat: 240664)
  - △ Dichlorodimethylsilane (Sigma, Cat: 440272)
- Acetone (Fisher, Cat: A18-1)

### Methods

1. Pour PDMS on SU-8 patterned master wafer; aluminum foil the bottom of the master for easy removal.
  - 1.1. Use 1:10 base-curing agent combination; 10-15 mL of total volume for 90 mm dishes, or 20-30 mL for 150 mm dishes work best for this application
2. Ink the stamps:
  - 2.1. Remove the elastomeric negative from the master and cut out the desired patterns from the whole stamp. For easy stamping, keep 0.5-1 cm pieces (0.25-1 cm<sup>2</sup>). *Optional:* Place the desired stamps on a plain glass microscope slide and plasma-clean for ~60 seconds just before use. This step is especially useful for stamping poly-L-lysine.
  - 2.2. Create a moisture-chamber by placing a 4"x4" Parafilm at the center of a 150mm tissue culture dish and adding pieces of wet KimWipes at the edges.
  - 2.3. Add desired number of protein solution droplets (~40 µL for a 1 cm<sup>2</sup> stamp) on the Parafilm and place the stamps upside-down into the droplets paying close attention to induce complete contact with the protein solution. Incubate at room temperature for 1 hour.

*Note:* If fluorescently labeled, keep in dark.
3. Prepare silanized glass coverslips:
  - 3.1. *Optional:* Immediately prior to stamping, plasma-clean the glass coverslips (that are already pre-cleaned with detergent and rinsed with diH<sub>2</sub>O and Ni dried) for ~60 seconds. This step is highly recommended for even and consistent silanization of surface.
  - 3.2. Place small amount of silanization solution into a glass dish (silane reacts with everything) under a chemical hood.

*Note:* Both benzene and silane are highly volatile and toxic; ensure proper safety.
  - 3.3. Add the coverslip into solution, swirl to ensure even coating, and incubate 20 seconds.
  - 3.4. Rapidly move the coverslip into acetone and wash for 10 seconds.
  - 3.5. Wash in running ddH<sub>2</sub>O for 30 seconds (as an option, use two consecutive diH<sub>2</sub>O baths prior to running water wash for complete removal of acetone and residual silane).

- 3.6. Blow dry coverslip with nitrogen.
4. Stamp the desired protein:
  - 4.1. Rinse the pattern-to-be-stamped with PBS, PBS-Tween, PBS and ddH<sub>2</sub>O (by dunking rapidly x30-50)
  - 4.2. Blow dry stamp with nitrogen without rewetting the surface.
  - 4.3. Put the glass coverslip on a large glass slide or culture dish; carefully grab the stamp and orient it at a negative angle (such that when it's upside down, forceps point down as well).
  - 4.4. Slowly place the stamp face down on the coverslip in one smooth motion. Use forceps to gently reinforce complete contact -- do not use excessive force since that may induce bending and over-contact between patterns.
  - 4.5. Gently place 12-50 g of weight (per cm<sup>2</sup>) on top of the stamp to establish firm contact between the PDMS stamp and the coverslip. A steel hexnut (3-4 g) is usually adequate for a standard 0.25 cm<sup>2</sup> stamp. If needed (i.e. poor transfer, pure edge contact, etc) use a steel bolt or thumbscrew (10-15 g).
  - 4.6. Wait for 30-120 seconds for adsorption of proteins on the coverslip. Transfer is likely to be an instant process; be consistent.
  - 4.7. Carefully remove the weight. *Optional:* Remove the stamp without any lateral motion; any motion may cause a smeared pattern or dislocation and destruction of transferred protein.
5. Coat the surface with Pluronic:
  - 5.1. Slowly add the Pluronic solution (without moving the stamp, if it is still on there) making sure it covers the entire coverslip. Incubate at room temperature for 1 hour. Remove the stamp after incubation.
  - 5.2. Carefully start washing with ddH<sub>2</sub>O while aspirating Pluronic; ensure that the surface is never exposed to air during this process. Use adequate water to remove all Pluronic; this may take 50-100 mL per coverslip.