



RNA Isolation with TRIzol

Prepared by: Alecia Muwonge, Nanditha Anandakrishnan

Version: 1.1

Date: 15-March-2022

Materials and equipment

- TRIzol (Life Technologies; Cat# 15596018)
- Cell culture dish containing cells
- 1X PBS (Boston BioProducts; Cat# BM-220), chilled on ice
- Chloroform (Sigma-Aldrich; Cat# C2432)
- 2-Propanol/Isopropyl Alcohol (Sigma Aldrich; Cat# I9516)
- Molecular Biology Grade Ethanol (Sigma- Aldrich; Cat# E7023) 70%, chilled on ice
- Nuclease Free Water (New England Biolabs Catalog # B1500S), chilled on ice
- Refrigerated microcentrifuge
- 2 mL microcentrifuge tubes
- DNA LoBind 1.5 mL Eppendorf tubes (Fisher Scientific; Cat# 022431021)
- Ice bucket with ice

Method

Cool PBS on ice. Place Culture dish containing cells on ice.



Discard medium. Quickly wash 2X with cold PBS.



Add 1 mL TRIzol into 60 mm dish. Incubate in hood for ≥ 5 minutes.



Lyse the cells directly in the culture dish by pipetting up and down several times



Transfer cells to 1.5 mL tubes



Centrifuge at 12,000g for 10 minutes at 4°C.



Transfer cleared supernatant to a new tube. Discard pellet.



Add 200 μ L Chloroform to the tube (200 μ L per 1 mL TRIzol).



Shake the tube vigorously by hand for 30 seconds.



Centrifuge at 12,000g for 10 minutes at 4°C.



Transfer the supernatant to a new tube. Avoid disturbing the interphase.



